

## Product Datasheet

# Porcine Epithelial Growth Factor ELISA Kit (orb1211920)

<b>Description</b>	Porcine Epithelial Growth Factor, EGF ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	200-20000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	EGF
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	133,506
<b>Uniprot ID</b>	<a href="#">Q00968</a>
<b>Solubility (25°C)</b>	24000ng/L
<b>Sensitivity</b>	98.23ng/L
<b>Entrez</b>	<a href="#">397083</a>
<b>Expiration Date</b>	Please enquire.