

Product Datasheet

Porcine Tumor Necrosis Factor A ELISA Kit (orb1211948)

Description	Porcine Tumor Necrosis Factor A, TNF-alpha ELISA Kit
Reactivity	Porcine
Range	10-4000ng/L
Tested Applications	ELISA
Target	TNFA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,254
Uniprot ID	P23563
Solubility (25°C)	4800ng/L
Sensitivity	5.06ng/L
Entrez	397086
Expiration Date	Please enquire.