

## Product Datasheet

# Porcine Insulin-like Growth Factors Binding Protein 3 ELISA Kit (orb1211949)

<b>Description</b>	Porcine Insulin-like Growth Factors Binding Protein 3, IGFBP3 ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	30-9000ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	IGFBP3
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	31,690
<b>Uniprot ID</b>	<a href="#">P16611</a>
<b>Solubility (25°C)</b>	9600ng/ml
<b>Sensitivity</b>	14.21ng/ml
<b>Entrez</b>	<a href="#">448812</a>
<b>Expiration Date</b>	Please enquire.