

## Product Datasheet

### Porcine Insulin ELISA Kit (orb1211953)

<b>Description</b>	Porcine Insulin, INS ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	3-900uIU/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	INS
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	11672
<b>Uniprot ID</b>	<a href="#">P01315</a>
<b>Solubility (25°C)</b>	960µIU/ml
<b>Sensitivity</b>	1.43µIU/ml
<b>Entrez</b>	<a href="#">397415</a>
<b>Expiration Date</b>	Please enquire.