

## Product Datasheet

# Porcine Endothelial Nitric Oxide Synthase ELISA Kit (orb1211957)

<b>Description</b>	Porcine Endothelial Nitric Oxide Synthase, ENOS ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	1-400IU/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	NOS3
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	133,406
<b>Uniprot ID</b>	<a href="#">Q28969</a>
<b>Solubility (25°C)</b>	480IU/mL
<b>Sensitivity</b>	0.55IU/mL
<b>Entrez</b>	<a href="#">397557</a>
<b>Expiration Date</b>	Please enquire.