

Product Datasheet

Porcine Cholecystokinin ELISA Kit (orb1211959)

Description	Porcine Cholecystokinin, CCK ELISA Kit
Reactivity	Porcine
Range	0.5-100ng/L
Tested Applications	ELISA
Target	CCK
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	12,526
Uniprot ID	P01356
Solubility (25°C)	128ng/L
Sensitivity	0.22ng/L
Entrez	397468
Expiration Date	Please enquire.