

## Product Datasheet

# Porcine Follicle-stimulating Hormone ELISA Kit (orb1211960)

<b>Description</b>	Porcine Follicle-stimulating Hormone, FSH ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	7-1500mIU/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	FSHR
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	78,173
<b>Uniprot ID</b>	<a href="#">P49059</a>
<b>Solubility (25°C)</b>	1600mIU/mL
<b>Sensitivity</b>	3.13mIU/mL
<b>Entrez</b>	<a href="#">397679</a>
<b>Expiration Date</b>	Please enquire.