

Product Datasheet

Porcine Interleukin 8 ELISA Kit (orb1211962)

Description	Porcine Interleukin 8, IL-8 ELISA Kit
Reactivity	Porcine
Range	1-380ng/L
Tested Applications	ELISA
Target	IL8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,633
Uniprot ID	P26894
Solubility (25°C)	400ng/L
Sensitivity	0.54ng/L
Entrez	396880
Expiration Date	Please enquire.