

Product Datasheet

Porcine Interleukin 4 ELISA Kit (orb1211964)

Description	Porcine Interleukin 4, IL-4 ELISA Kit
Reactivity	Porcine
Range	1-380ng/L
Tested Applications	ELISA
Target	IL4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	15,023
Uniprot ID	Q04745
Solubility (25°C)	400ng/L
Sensitivity	0.51ng/L
Entrez	397225
Expiration Date	Please enquire.