

Product Datasheet

Porcine Interleukin 1A ELISA Kit (orb1211967)

Description	Porcine Interleukin 1A, IL-1A ELISA Kit
Reactivity	Porcine
Range	0.5-300ng/L
Tested Applications	ELISA
Target	IL1A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,789
Uniprot ID	P18430
Solubility (25°C)	320ng/L
Sensitivity	0.33ng/L
Entrez	397094
Expiration Date	Please enquire.