

Product Datasheet

Porcine Interleukin 18 ELISA Kit (orb1211968)

Description	Porcine Interleukin 18, IL-18 ELISA Kit
Reactivity	Porcine
Range	20-4500ng/L
Tested Applications	ELISA
Target	IL18
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,026
Uniprot ID	O19073
Solubility (25°C)	4800ng/L
Sensitivity	9.66ng/L
Entrez	397057
Expiration Date	Please enquire.