

Product Datasheet

Porcine Interleukin 13 ELISA Kit (orb1211969)

Description	Porcine Interleukin 13, IL-13 ELISA Kit
Reactivity	Porcine
Range	5-600ng/L
Tested Applications	ELISA
Target	IL13
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,452
Uniprot ID	Q95J68
Solubility (25°C)	640ng/L
Sensitivity	2.31ng/L
Entrez	396721
Expiration Date	Please enquire.