

## Product Datasheet

### Porcine Interleukin 10 ELISA Kit (orb1211970)

<b>Description</b>	Porcine Interleukin 10, IL-10 ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	3-600ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	IL10
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	19,934
<b>Uniprot ID</b>	<a href="#">Q29055</a>
<b>Solubility (25°C)</b>	640ng/L
<b>Sensitivity</b>	1.46ng/L
<b>Entrez</b>	<a href="#">397106</a>
<b>Expiration Date</b>	Please enquire.