

Product Datasheet

Porcine E-Selectin ELISA Kit (orb1211972)

| | |
|----------------------------|--|
| Description | Porcine E-Selectin ELISA Kit |
| Reactivity | Porcine |
| Range | 1-400ng/mL |
| Tested Applications | ELISA |
| Target | SELE |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 52,567 |
| Uniprot ID | P98110 |
| Solubility (25°C) | 480ng/mL |
| Sensitivity | 0.56ng/mL |
| Entrez | 397508 |
| Expiration Date | Please enquire. |