

Product Datasheet

Mouse Nad-dependent Protein Deacetylase Sirtuin-6 ELISA Kit (orb1212013)

Description	Mouse Nad-dependent Protein Deacetylase Sirtuin-6, SIRT6 ELISA Kit
Reactivity	Mouse
Range	7-1500ng/L
Tested Applications	ELISA
Target	SIRT6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,920
Uniprot ID	P59941
Solubility (25°C)	1600ng/L
Sensitivity	3.68ng/L
Entrez	50721
Expiration Date	Please enquire.