

Product Datasheet

Mouse Synaptophysin ELISA Kit (orb1212016)

Description	Mouse Synaptophysin, SYP ELISA Kit
Reactivity	Mouse
Range	0.5-300ng/ml
Tested Applications	ELISA
Target	SYP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	34,025
Uniprot ID	Q62277
Solubility (25°C)	320ng/ml
Sensitivity	0.3ng/ml
Entrez	20977
Expiration Date	Please enquire.