

Product Datasheet

Mouse Tumor Necrosis Factor Ligand Superfamily Member 13 ELISA Kit (orb1212029)

Description	Mouse Tumor Necrosis Factor Ligand Superfamily Member 13, TNFSF13 ELISA Kit
Reactivity	Mouse
Range	50-10000ng/L
Tested Applications	ELISA
Target	TNFSF13
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,889
Uniprot ID	Q9D777
Solubility (25°C)	12800ng/L
Sensitivity	24.83ng/L
Entrez	69583
Expiration Date	Please enquire.