

Product Datasheet

Mouse Sortilin ELISA Kit (orb1212047)

Description	Mouse Sortilin, SORT1 ELISA Kit
Reactivity	Mouse
Range	3-900ng/L
Tested Applications	ELISA
Target	SORT1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	91,200
Uniprot ID	Q6PHU5
Solubility (25°C)	960ng/L
Sensitivity	1.46ng/L
Entrez	20661
Expiration Date	Please enquire.