

Product Datasheet

Mouse Semaphorin-4D ELISA Kit (orb1212049)

Description	Mouse Semaphorin-4D, SEMA4D ELISA Kit
Reactivity	Mouse
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	SEMA4D
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	95,640
Uniprot ID	O09126
Solubility (25°C)	48ng/mL
Sensitivity	0.056ng/mL
Entrez	20354
Expiration Date	Please enquire.