

Product Datasheet

Mouse R-spondin-1 ELISA Kit (orb1212059)

Description	Mouse R-spondin-1, RSPO1 ELISA Kit
Reactivity	Mouse
Range	2-700ng/L
Tested Applications	ELISA
Target	RSPO1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,332
Uniprot ID	Q9Z132
Solubility (25°C)	800ng/L
Sensitivity	0.98ng/L
Entrez	192199
Expiration Date	Please enquire.