

## Product Datasheet

### Mouse Prokineticin-2 ELISA Kit (orb1212075)

<b>Description</b>	Mouse Prokineticin-2, PROK2 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	7-1500ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	PROK2
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	14,185
<b>Uniprot ID</b>	<a href="#">Q9QXU7</a>
<b>Solubility (25°C)</b>	1600ng/L
<b>Sensitivity</b>	3.43ng/L
<b>Entrez</b>	<a href="#">50501</a>
<b>Expiration Date</b>	Please enquire.