

Product Datasheet

Mouse Programmed Cell Death 1 Ligand 1 ELISA Kit (orb1212076)

Description	Mouse Programmed Cell Death 1 Ligand 1, CD274 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	CD274
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	32,780
Uniprot ID	Q9EP73
Solubility (25°C)	4800ng/L
Sensitivity	9.21ng/L
Entrez	60533
Expiration Date	Please enquire.