

Product Datasheet

Mouse Peptide Yy ELISA Kit (orb1212087)

Description	Mouse Peptide Yy, PYY ELISA Kit
Reactivity	Mouse
Range	50-10000ng/L
Tested Applications	ELISA
Target	PYY
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,064
Uniprot ID	Q9EPS2
Solubility (25°C)	12800ng/L
Sensitivity	23.35ng/L
Entrez	217212
Expiration Date	Please enquire.