

Product Datasheet

Mouse Pentraxin-related Protein Ptx3 ELISA Kit (orb1212088)

Description	Mouse Pentraxin-related Protein Ptx3, PTX3 ELISA Kit
Reactivity	Mouse
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	PTX3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	41,812
Uniprot ID	P48759
Solubility (25°C)	48ng/mL
Sensitivity	0.052ng/mL
Entrez	19288
Expiration Date	Please enquire.