

Product Datasheet

Mouse Neuropeptide S ELISA Kit (orb1212104)

Description	Mouse Neuropeptide S, NPS ELISA Kit
Reactivity	Mouse
Range	2-600ng/L
Tested Applications	ELISA
Target	NPS
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	9,861
Uniprot ID	P0C0P8
Solubility (25°C)	640ng/L
Sensitivity	0.97ng/L
Entrez	100043254
Expiration Date	Please enquire.