

## Product Datasheet

# Mouse Neuropeptide Ff Receptor 2 ELISA Kit (orb1212105)

<b>Description</b>	Mouse Neuropeptide Ff Receptor 2, NPFFR2 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	20-4500ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	NPFFR2
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	47,449
<b>Uniprot ID</b>	<a href="#">Q924H0</a>
<b>Solubility (25°C)</b>	4800ng/L
<b>Sensitivity</b>	9.81ng/L
<b>Entrez</b>	<a href="#">104443</a>
<b>Expiration Date</b>	Please enquire.