

Product Datasheet

Mouse Muscarinic Acetylcholine Receptor M1 ELISA Kit (orb1212117)

Description	Mouse Muscarinic Acetylcholine Receptor M1, CHRM1 ELISA Kit
Reactivity	Mouse
Range	10-2000ng/L
Tested Applications	ELISA
Target	CHRM1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,379
Uniprot ID	P12657
Solubility (25°C)	2400ng/L
Sensitivity	5.09ng/L
Entrez	12669
Expiration Date	Please enquire.