

Product Datasheet

Mouse Metastasis-Suppressor Kiss-1 ELISA Kit (orb1212128)

Description	Mouse Metastasis-Suppressor Kiss-1, KISS1 ELISA Kit
Reactivity	Mouse
Range	7-1500ng/L
Tested Applications	ELISA
Target	KISS1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,117
Uniprot ID	Q6Y4S4
Solubility (25°C)	1600ng/L
Sensitivity	3.45ng/L
Entrez	280287
Expiration Date	Please enquire.