

## Product Datasheet

### Mouse Metallothionein-2 ELISA Kit (orb1212129)

<b>Description</b>	Mouse Metallothionein-2, MT-2 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	10-2000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	MT2
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	6,115
<b>Uniprot ID</b>	<a href="#">P02798</a>
<b>Solubility (25°C)</b>	2400ng/L
<b>Sensitivity</b>	5.54ng/L
<b>Entrez</b>	<a href="#">17750</a>
<b>Expiration Date</b>	Please enquire.