

Product Datasheet

Mouse Lipopolysaccharide-binding Protein ELISA Kit (orb1212149)

Description	Mouse Lipopolysaccharide-binding Protein, LBP ELISA Kit
Reactivity	Mouse
Range	0.3-90ug/ml
Tested Applications	ELISA
Target	LBP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,055
Uniprot ID	Q61805
Solubility (25°C)	96ug/ml
Sensitivity	0.12ug/ml
Entrez	16803
Expiration Date	Please enquire.