

Product Datasheet

Mouse Interleukin-19 ELISA Kit (orb1212169)

Description	Mouse Interleukin-19, IL19 ELISA Kit
Reactivity	Mouse
Range	7-1500ng/L
Tested Applications	ELISA
Target	IL19
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,288
Uniprot ID	Q8CJ70
Solubility (25°C)	1600ng/L
Sensitivity	3.16ng/L
Entrez	329244
Expiration Date	Please enquire.