

Product Datasheet

Mouse Heat Shock Cognate 71 kDa Protein ELISA Kit (orb1212205)

Description	Mouse Heat Shock Cognate 71 kDa Protein, HSPA8 ELISA Kit
Reactivity	Mouse
Range	5-1000ng/L
Tested Applications	ELISA
Target	HSPA8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	70,871
Uniprot ID	P63017
Solubility (25°C)	1280ng/L
Sensitivity	2.45ng/L
Entrez	15481
Expiration Date	Please enquire.