

## Product Datasheet

### Mouse Glia-derived Nexin ELISA Kit (orb1212218)

<b>Description</b>	Mouse Glia-derived Nexin, SERPINE2 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	1-400ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	SERPINE2
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	44,207
<b>Uniprot ID</b>	<a href="#">Q07235</a>
<b>Solubility (25°C)</b>	480ng/L
<b>Sensitivity</b>	0.56ng/L
<b>Entrez</b>	<a href="#">20720</a>
<b>Expiration Date</b>	Please enquire.