

Product Datasheet

Mouse Ephrin-A1 ELISA Kit (orb1212239)

Description	Mouse Ephrin-A1, EFNA1 ELISA Kit
Reactivity	Mouse
Range	50-10000ng/L
Tested Applications	ELISA
Target	EFNA1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,802
Uniprot ID	P52793
Solubility (25°C)	12800ng/L
Sensitivity	23.43ng/L
Entrez	13636
Expiration Date	Please enquire.