

Product Datasheet

Mouse EGF-like Module-containing Mucin-like Hormone Receptor-like 1 ELISA Kit (orb1212245)

Description	Mouse EGF-like Module-containing Mucin-like Hormone Receptor-like 1, EMR1 ELISA Kit
Reactivity	Mouse
Range	3-900ng/L
Tested Applications	ELISA
Target	ADGRE1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	102,130
Uniprot ID	Q61549
Solubility (25°C)	960ng/L
Sensitivity	1.65ng/L
Entrez	13733
Expiration Date	Please enquire.