

Product Datasheet

Mouse EGF-containing Fibulin-like Extracellular Matrix Protein 1 ELISA Kit (orb1212246)

Description	Mouse EGF-containing Fibulin-like Extracellular Matrix Protein 1, EFEMP1 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	EFEMP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,953
Uniprot ID	Q8BPB5
Solubility (25°C)	4800ng/L
Sensitivity	9.95ng/L
Entrez	216616
Expiration Date	Please enquire.