

## Product Datasheet

# Mouse Disintegrin and Metalloproteinase Domain-containing Protein 12 ELISA Kit (orb1212255)

<b>Description</b>	Mouse Disintegrin and Metalloproteinase Domain-containing Protein 12, ADAM12 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.05-20ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	ADAM12
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	98,504
<b>Uniprot ID</b>	<a href="#">Q61824</a>
<b>Solubility (25°C)</b>	24ng/mL
<b>Sensitivity</b>	0.024ng/mL
<b>Entrez</b>	<a href="#">11489</a>
<b>Expiration Date</b>	Please enquire.