

## Product Datasheet

# Mouse Dickkopf-related Protein 3 ELISA Kit (orb1212265)

<b>Description</b>	Mouse Dickkopf-related Protein 3, DKK3 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	5-1000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	DKK3
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	38,388
<b>Uniprot ID</b>	<b>Q9QUN9</b>
<b>Solubility (25°C)</b>	1280ng/L
<b>Sensitivity</b>	2.35ng/L
<b>Entrez</b>	<b>50781</b>
<b>Expiration Date</b>	Please enquire.