

Product Datasheet

Mouse Dickkopf-related Protein 1 ELISA Kit (orb1212266)

Description	Mouse Dickkopf-related Protein 1, DKK1 ELISA Kit
Reactivity	Mouse
Range	30-9000ng/L
Tested Applications	ELISA
Target	DKK1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,298
Uniprot ID	O54908
Solubility (25°C)	9600ng/L
Sensitivity	14.12ng/L
Entrez	13380
Expiration Date	Please enquire.