

## Product Datasheet

# Mouse Nt-3 Growth Factors Receptor ELISA Kit (orb1212312)

<b>Description</b>	Mouse Nt-3 Growth Factors Receptor, NTRK2 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	20-4500ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	NTRK2
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	92,133
<b>Uniprot ID</b>	<b>P15209</b>
<b>Solubility (25°C)</b>	4800ng/ml
<b>Sensitivity</b>	9.23ng/ml
<b>Entrez</b>	<b>18212</b>
<b>Expiration Date</b>	Please enquire.