

Product Datasheet

Mouse C-C Motif Chemokine 2 ELISA Kit (orb1212346)

| | |
|----------------------------|---|
| Description | Mouse C-C Motif Chemokine 2, CCL2 ELISA Kit |
| Reactivity | Mouse |
| Range | 2-600pg/ml |
| Tested Applications | ELISA |
| Target | CCL2 |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 16,326 |
| Uniprot ID | P10148 |
| Solubility (25°C) | 800pg/ml |
| Sensitivity | 0.93pg/ml |
| Entrez | 20296 |
| Expiration Date | Please enquire. |