

Product Datasheet

Mouse Toll-like Receptor 4 ELISA Kit (orb1212351)

Description	Mouse Toll-like Receptor 4, TLR-4 ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	TLR4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	95,519
Uniprot ID	Q9QUK6
Solubility (25°C)	24ng/ml
Sensitivity	0.023ng/ml
Entrez	21898
Expiration Date	Please enquire.