

Product Datasheet

Mouse Caspase 3 ELISA Kit (orb1212380)

Description	Mouse Caspase 3, CASP3 ELISA Kit
Reactivity	Mouse
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	CASP3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	31475
Uniprot ID	P70677
Solubility (25°C)	12.8ng/ml
Sensitivity	0.023ng/ml
Entrez	12367
Expiration Date	Please enquire.