

Product Datasheet

Mouse TGF-beta Receptor Type-1 ELISA Kit (orb1212396)

Description	Mouse TGF-beta Receptor Type-1, TGFBR1 ELISA Kit
Reactivity	Mouse
Range	3-900ng/L
Tested Applications	ELISA
Target	TGFBR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	56,179
Uniprot ID	Q64729
Solubility (25°C)	960ng/L
Sensitivity	1.53ng/L
Entrez	21812
Expiration Date	Please enquire.