

Product Datasheet

Mouse Signal Peptide ELISA Kit (orb1212439)

Description	Mouse Signal Peptide, Cub and EGF-like Domain-containing Protein 1, SCUBE1 ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	SCUBE1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	111,606
Uniprot ID	Q6NZL8
Solubility (25°C)	24ng/mL
Sensitivity	0.027ng/mL
Entrez	64706
Expiration Date	Please enquire.