

Product Datasheet

Mouse Interleukin 17A ELISA Kit (orb1212452)

Description	Mouse Interleukin 17A, IL-17A ELISA Kit
Reactivity	Mouse
Range	2-600ng/L
Tested Applications	ELISA
Target	IL17A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,490
Uniprot ID	Q62386
Solubility (25°C)	800ng/L
Sensitivity	0.97ng/L
Entrez	16171
Expiration Date	Please enquire.