

Product Datasheet

Mouse Secretin ELISA Kit (orb1212453)

Description	Mouse Secretin, SCT ELISA Kit
Reactivity	Mouse
Range	7-1500ng/L
Tested Applications	ELISA
Target	SCT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,914
Uniprot ID	Q08535
Solubility (25°C)	1600ng/L
Sensitivity	3.43ng/L
Entrez	20287
Expiration Date	Please enquire.