

Product Datasheet

Mouse Ghrelin ELISA Kit (orb1212460)

Description	Mouse Ghrelin, GHRL ELISA Kit
Reactivity	Mouse
Range	1-400ng/L
Tested Applications	ELISA
Target	GHRL
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	13,207
Uniprot ID	Q9EQX0
Solubility (25°C)	480ng/L
Sensitivity	0.58ng/L
Entrez	58991
Expiration Date	Please enquire.