

## Product Datasheet

### Mouse Apelin ELISA Kit (orb1212473)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Mouse Apelin, AP ELISA Kit             |
| <b>Reactivity</b>          | Mouse                                  |
| <b>Range</b>               | 5-2000pg/ml                            |
| <b>Tested Applications</b> | ELISA                                  |
| <b>Target</b>              | APLN                                   |
| <b>Storage</b>             | 2-8°C                                  |
| <b>Note</b>                | For research use only                  |
| <b>Application notes</b>   | Sandwich                               |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates |
| <b>MW</b>                  | 8,658                                  |
| <b>Uniprot ID</b>          | <a href="#">Q9R0R4</a>                 |
| <b>Solubility (25°C)</b>   | 2400pg/ml                              |
| <b>Sensitivity</b>         | 3.11pg/ml                              |
| <b>Entrez</b>              | <a href="#">30878</a>                  |
| <b>Expiration Date</b>     | Please enquire.                        |