

Product Datasheet

Mouse Cyclooxygenase-1 ELISA Kit (orb1212476)

Description	Mouse Cyclooxygenase-1, COX-1 ELISA Kit
Reactivity	Mouse
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	PTGS1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	69,042
Uniprot ID	P22437
Solubility (25°C)	80ng/ml
Sensitivity	0.13ng/ml
Entrez	19224
Expiration Date	Please enquire.